

SM-8 Electronic Load Control

STERILIZATION RECORD KEEPING LABELING SYSTEM

Stand Alone Load Control System with Barcode Support


Validated for ALL STERRAD Systems



- BARCODE INFORMATION COMPLIES WITH **REAL TIME TRACEABILITY** HOSPITAL REQUIREMENTS
- STAND-ALONE OR PC COMPATIBLE
- STEAM, GAS & STERRAD VALIDATED LABELS
- AUTOMATIC TRAY OR BUNDLE COUNTER
- INTERNAL MEMORY FOR 3,000 ITEMS
- TECHNICIAN IDENTIFICATION AREA
- HOSPITAL AND DEPARTMENT IDENTIFICATION
- PRE-PROGRAMMED - READY TO BE USED OUT OF THE BOX




CATALOG NUMBER: L51-A55EY
SYNTHETIC



CATALOG NUMBER: L51-A55EB
SYNTHETIC



CATALOG NUMBER: L51-A55W
SYNTHETIC



CATALOG NUMBER: L51-P55W
PAPER

OPTIONAL PRINTS - COLOR LABELS AVAILABLE

SM-8 Sterilization Label Coder Kit

The SM8 Load Control Labeling System offers the latest technology for your sterilization record keeping application. It brings all your load control label information plus barcode functionality and full tray or bundle item description.

Barcodes now can be scanned when the item is deployed to its destination minimizing mistakes in the information transfer by hand. The system offers enough memory to save all your inventory items and the tray or bundle information is readily available after scanning your tray library book or your tray ID barcode. The system asks for your new load control number, the technician ID and how many labels to print.

Custom layouts to maximize your labeling needs are available with paper and synthetic material options to support all sterilization methods including your STERRAD system.

For more information, brochure or sample labels printed from this new revolutionary system

contact us at **1.800.679.7836**
fax: 813-972-0578 - www.S1Medical.com

SM8 LOAD CONTROL SYSTEM



CATALOG NO.	DESCRIPTION
Z70-S8MKT	SM-8 STAND ALONE FULL KIT -PRINTER, KEYBOARD, CABLE, SOFTWARE -20 DOWNLOADS
L51-P55 SERIES LABELS	PAPER - LABEL 1.75" x 1" 8 ROLLS / BOX - COLORS: WHITE, YELLOW, PINK & BLUE RIBBON: *R10-ZP17 - 12 ROLLS PER BOX
L51-A55 SERIES LABELS	SYNTHETIC - LABEL 1.75" x 1" 8 ROLLS / BOX - COLORS: WHITE, YELLOW, PINK & BLUE RIBBON: *R10-ZS17 - 12 ROLLS PER BOX
A70-CSKIT	SM-1 CODER KIT - 100+ PRE-LOADED LABELS
M91-DS4800c	CORDED SCANNER WITH CABLE ADAPTER