




# Central Sterile Labeling & Tracking Systems

**SteriTag™ System**  
Color Tracking Sterilization Tags







## Create Your Custom SteriTags™ Safely Customize Your Trays With Custom Colorful Tags!

- CUSTOMIZE YOUR TAG COLOR and INFORMATION
- FULLY COMPATIBLE WITH ALL TRAY TRACKING SYSTEMS
- FLAT or POCKET COLOR TAGS AVAILABLE
- COLOR IDENTIFY ALL YOUR TRAYS and ITEMS
- ASP FULL VALIDATION (for all 5 STERRAD Systems)
- STEAM, GAS & STERRAD COMPATIBLE
- BRILLIANT COLORS - WRITE FRIENDLY TAGS
- SYNTHETIC STERILIZATION LASER SHEETS- AVAILABLE
- FAST DELIVERY ON ALL CUSTOM TAGS



**SteriPocket™ System**  
(pocket color tags)  
shown with synthetic sheet  
Cat. #: L51-T005SY

**SteriTag™ System**  
(flat color tags)  
Cat. #: L51-T005 Stock



Our SteriTag™ System solution can be customized to meet all of your particular coding demands. Two formats are available to begin your new custom tag system:

- SteriTags** (flat surface color tags)
- SteriPockets** (pocket cuts color tags)

Our sterilization tags allow you to handwrite directly on to the tag or you can apply your information label to further enhance your current system. SteriTags™ are manufactured with our exclusive synthetic material that have been fully tested for all sterilization processes such as Steam, ETO and the complete STERRAD System Line (ASP SSG Guide Validated) including all washing procedures.



## Central Sterile Labeling & Tracking

# SteriTag™ System

Tracking Sterilization Tags

The Cost Effective Solution for Color Coding Your Instruments, Vendor and Repair Trays and more...




- INDIVIDUALLY IDENTIFY ALL YOUR TRAYS and ITEMS
- WRITE FRIENDLY SURFACE
- ASP VALIDATED (Advanced Sterilization Products)
- STEAM, GAS & STERRAD COMPATIBLE
- SIX COLORS IN STOCK
- CUSTOM PRINT AND COLORS AVAILABLE
- FULLY COMPATIBLE WITH ALL TRACKING SYSTEM APPLICATIONS
- CONVENIENT 5 COLOR KIT-100 PER COLOR

System One Medical is proud to announce our New SteriTag™ Tracking System. This is the most cost effective solution to identifying all your trays fast and efficiently. Our SteriTag™ Tracking System solution offers several colors to choose from and they are convenient to write directly on the tag or place your barcode label directly on it for complete electronic tray identification.



Our SteriTags™ have been specifically designed to improve item management, tray identification and reduce costly mistakes. This new system dramatically enhances your tray identification and tray maintenance schedule proceedings. SteriTags™ are manufactured with our exclusive synthetic material that have been fully tested for all sterilization procedures such as Steam, ETO and STERRAD (ASP Validated) including all washing actions. These versatile tags are fully compatible with all bar code labeling system including Applied Logic / Abacus, Microsystems, Cardinal, Censitrac, T-DOC, Jarit, ABGI, Rosebud, etc.

Without delay, you can start your secure and cost effective item management solution today. Simply write your information directly on the tag or place your tracking label over your SteriTag™ color of choice over and over again! Custom prints, additional colors and designs are available.



TAG SIZE: 8" x 2.25"		TAGS PER BOX: 500
L51-T005W	SteriTAG™ - WHITE	
L51-T005Y	SteriTAG™ - FLO-YELL	
L51-T005B	SteriTAG™ - FLO-BLUE	
L51-T005P	SteriTAG™ - FLO-PINK	
L51-T005R	SteriTAG™ - FLO-RED	
L51-T005G	SteriTAG™ - FLO-GREEN	

For more information, data sheets or samples call us or visit us at  
**1.800.679.7836**    **WWW.S1MEDICAL.COM**

**System One** MEDICAL  
 BARCODE AND LABEL SYSTEMS SPECIALISTS